

# Rabies Serology



Greendale Veterinary  
Diagnostics

Delivering Results



**BioBest**

*Veterinary Serology at its Best*

- Use one submission form per animal.
- Send a minimum of 1ml serum (preferable) or 2ml clotted blood.
- Label sample clearly with the animal's microchip number.
- For the purpose of the UK Pet Travel Scheme, a test titre result of 0.5IU/ml or above indicates that the dog or cat has an acceptable rabies antibody level.
- Sample will be tested by BioBest.

Greendale Veterinary Diagnostics Ltd,  
Lansbury Estate,  
Knaphill,  
Woking,  
Surrey GU21 2EW

Greendale Lab No:

## Owner's Details

Name:

Address:

## Submitting Veterinary Surgeon's Details

Practice Name and Address:

Signature of submitting veterinary surgeon:

Tel No:

Name (caps):

Fax No:

Date:



### Rabies Extra Service

Tick box to include a diagnostic profile for only £4.00 extra  
(TP, Alb, Glob, ALT, ALP, Urea, Creatinine, Total Bilirubin (if icteric))

## Animal's Details

Microchip Number:

Animal Name:

Date of Birth:

Species:

Date of Sampling and Microchip reading:

Rabies Vaccination Details:

Date	Vaccine	Batch No.

## For BioBest Use Only

Form V1.1

Date Received:

QC:

F:

I:

BioBest Ref: